



**Finnish Institute of
Occupational Health**



BRAIN•WORK
RESEARCH CENTRE

MIDAS

Framework for Distributed Online Data Stream Analysis

Andreas Henelius andreas.henelius@ttl.fi
Jari Torniainen jari.torniainen@ttl.fi

Finnish Institute of Occupational Health
Brain Work Research Centre
Helsinki, Finland

Symbiotic2015, Berlin, 07.10.2015

OUTLINE

Part I (1h)

- Background
- Architecture (high-level)
- The MIDAS node
- Configuration and usage
- Building analysis nodes
- Installation
- Resources

Part II (2h)

- Practical Exercises

M I D A S


MODULAR INTEGRATED DISTRIBUTED ANALYSIS SYSTEM

WHAT IS MIDAS?

A system for online analysis of streaming signals and allows easy integration of such into machine learning frameworks.

WHY ONLINE ANALYSIS?

- Cognitive State Determination
- Symbiotic Interaction
- IoT



CHALLENGES


- Sensors and protocols
- Data streams and formats
- Signal fusion and computational load

WHAT IS MIDAS?

Modular Integrated Distributed Analysis System




SIMILAR SOLUTIONS




OpenViBE

Renard, Yann, et al. "OpenViBE: an open-source software platform to design, test, and use brain-computer interfaces in real and virtual environments" Presence: teleoperators and virtual environments, 2010




SSI

Wagner, Johannes, et al. "The Social Signal Interpretation Framework (SSI) for Real Time Signal Processing and Recognition" INTERSPEECH, 2011



BCILAB

C Kothe, S Makeig. "BCILAB: A platform for brain-computer interface development" Journal of Neural Engineering, 2013



BCI2000

G. Schalk, D.J. McFarland, T. Hinterberger, N. Birbaumer, and J.R. Wolpaw: BCI2000: A General-Purpose Brain-Computer Interface (BCI) System, IEEE Trans Biomed Eng, 51(6), June 2004.

HOW IS MIDAS DIFFERENT?

Easily integrate and combine different sensors into online analyses of streaming data.


- **Slim:** minimal clutter
- **Modular:** add functionality as needed
- **Distributed:** Runs everywhere

- Written in Python
- Open Source (MIT)


MIDAS ARCHITECTURE




MIDAS ARCHITECTURE




MIDAS ARCHITECTURE



COMMUNICATION PROTOCOLS IN MIDAS




THE NODE




STREAMERS

- Read sensor data and retransmit raw data using the **lab streaming layer**
- Some devices have built-in LSL support (Neuroelectrics Enobio)
- For some devices you can find LSL drivers from **<http://code.google.com/p/labstreaminglayer/>**
- If no streamer exist, write it yourself! (not difficult)

BUILDING NODES




DATA CONTAINERS



- Receiver automatically stores all incoming data into primary data buffers
- Secondary data buffers can be used to store values calculated from the primary data
- Both primary and secondary data can be requested through the dispatcher
- Both primary and secondary data can be used as an input for the analysis functions

ANALYSIS FUNCTIONS



Nodes in the MIDAS framework have two types of analysis functions: **metrics** and **processes**

METRIC FUNCTION

- $y = f(x)$
- Calculates 'a metric' using data as input
- Input data is specified when the request is made by the client
- Return value can be a scalar or a vector/array
- Can be a method of the Node-class or an external function


METRIC FUNCTION

- The first positional argument of the metric function always contains the data
- The incoming data is wrapped inside a dict:

| Keys | Values | | | |
|--------|------------|------------|-----|------------|
| 'data' | samples | samples | ... | samples |
| 'time' | timestamps | timestamps | ... | timestamps |
| | channel 0 | channel 1 | | channel N |

PROCESS

- Automatically calculate values from the primary data at set intervals
- Runs in a separate process
- Pushes calculated values into secondary data buffer
- Must be a method of the Node-class



CODE EXAMPLE


Example from the MIDAS repository.

CONFIGURATION AND USE

1. Configure options in the INI-file
2. Ensure that the data stream is visible
3. Start the node(s)
4. Start the dispatcher
5. Access nodes through the MIDAS API (= client)

CONFIGURATION FILE

CONFIG.INI



CONFIGURATION FILE

Dispatcher

```
# Dispatcher-section
[dispatcher]
port                = 8080
ip                  = localhost
n_threads           = 5
run_pubsub_proxy    = False
proxy_port_in       = 5999
proxy_port_out      = 6000
```

CONFIGURATION FILE

An ECG node

```
# Real-time ECG node
[ecg_rt]
nodename           = ecgnode
nodetype           = ECG
nodeid             = 05
nodedesc           = ECG Processing node
primary_node       = True
port_frontend      = 5054
port_backend       = 5055
port_publisher     = 5056
run_publisher      = False
n_workers          = 2
n_channels         = 1
channel_names      = ch0
channel_descriptions = ECG channel
sampling_rate      = 100
buffer_size_s      = 600
lsl_stream_name    = faros_ecg
```

STARTING MIDAS

The components of MIDAS are started from the command line.

Dispatcher

```
midas-dispatcher config.ini dispatcher
```

A node

```
python3 ecg_node.py config.ini ecg  
./ecg_node.py config.ini ecg      [if first line is #!/usr/bin/env pytho  
n3]
```

MIDAS STATUS

Check the status of the MIDAS network

<http://127.0.0.1:8080>

<http://127.0.0.1:8080/status/nodes>

<http://127.0.0.1:8080/status/metrics>

QUERYING DATA

Python

```
import requests  
r = requests.get("http://localhost:8080/example_node/  
metric/{\"type\":\"metric_a\", \"time_window\":[5]}")
```

MATLAB

```
r = urlread('http://localhost:8080/example_node/  
metric/{\"type\":\"metric_a\", \"time_window\":[5]}');
```

MIDAS also supports **JSONP** for use with JavaScript.

INSTALLATION

Remember to use Python 3!

```
pip install git+https://github.com/bwrc/midas
```


RESOURCES

Source and Wiki

github.com/bwrc/midas

The API

github.com/bwrc/midas/wiki/API

FUTURE

Mobile Applications: Move experiments out of the lab

Internet of Things: Utilize wearable health technology

CONTRIBUTE!

Together we can make MIDAS what it needs to be!

Q & A